

Power voltage 230 V, 50/60 Hz

Batteries Batteries 2 x 12 V / 28 Ah
Input power – mains 380 W
External source 195 W

Temperature adjustment range 28.0 – 38.0 °C continuously
O₂ concentration adjustment range 21 – 50 %
Inside noise maximum level max. 60 dB

Dimensions

Height 681 mm
Length 867 mm
Width 591 mm

The following conditions are signalized by alarms

- power failure
- high/low air temperature
- ventilator failure
- probe failure
- overheating
- heater disconnection
- power drop

Operation time with fully charged batteries min. 4 hrs

Optional accessories

Neotic Extero type carriage for ambulances
PNT – 401 type carriage for transport within hospitals
Rail with a shelf and infusion holder
LED examination light
Fixing belts
Cotton blanket for covering the patient space