

TSE

Smart technologies



BABY BLUE LIGHT FTL – 600

Phototherapy lamp

Baby Blue Light FTL – 600 Phototherapy lamp

The phototherapy lamp Baby Blue Light FTL – 600 represents the most common method of adequate treatment of neonatal hyperbilirubinemia. It brings an economical solution of adequate treatment along with a perfect price/performance ratio.

The simple suspension of the lamp on a height-adjustable stand enables both tilting up to 90° in one direction, and effective irradiation of the patient's entire body on an infant bed or a warming bed. Anti slip pads on the lamp surface ensure safe location on an incubator cover after removing the lamp from the stand.

A counter with a digital display records the number of operation hours and thus ensures safe and efficient operation.

Power voltage	230 V, 50 Hz
Max. Input power of the unit	170 VA
Input power of light sources	2 x 20 W
Light source	discharge tube PHILIPS (blue)
Wave length	400 – 550 nm
Light source life	approx. 2 000 hours
Dimensions	
Height	115 mm
Length	820 mm
Width	210 mm
Weight	
	4.5 kg
Optional accessories	
Stand FTL – 600	adjustable height 1025 – 1475 mm



TSE spol. s r. o.
Mánesova 74, 371 52 České Budějovice, Czech Republic
T: + 420 386 721 149, F: + 420 386 721 102
medical@tse.cz, www.tse.cz

