

**Power voltage**

9V ss

**Number of measurements  
with fully loaded battery**

720 measurements

**Dimensions without sensor:**

height

27 mm

length

110,5 mm

width

66,2 mm

**Measured light wavelength**

400 – 550 nm

**Measurement accuracy**

$\pm 10\% \pm 2$ digit.

**Measurement range**

Peak value (PEAK): 10 – 115  $\mu\text{W}/\text{cm}^2/\text{nm}$

Mean value (MEAN): 3 – 46  $\mu\text{W}/\text{cm}^2/\text{nm}$

**Calibration validity**

1 year